JUZO COMPRESSION WRAP

Care & Application for Juzo Compression Wrap
Foot Segment

REVERSIBLE

beige or black garment

Freedom in Motion

Thank you for choosing Juzo, an international leader in medical compression garments since 1912.

Manufactured with care and detail, your Juzo Compression Wrap has several features that provide therapy, comfort, versatility and superior fit, including:

- Working pressure - Short stretch for a high working pressure.
- Gap resistant - Special notches and minimal overlying straps resist gaping and bulk.
- EZ pull tab - Aid for hassle free application.
- Reversible - Worn in either beige or black

Included In This Package:

- Juzo Compression Wrap
- Juzo Compression Cotton Stockinette
- Juzo Hook Fasteners

Care and Maintenance

- Remove all the half circle hook fasteners from the straps
- Machine wash warm
- A garment bag is strongly recommended
- Tumble dry low
- Do not bleach
- DO NOT USE: bleach, chlorine, fabric softeners, or other laundry additives
- Store garments in a dry place at room temperature

Return Policy

Should it be necessary to return your garment for evaluation, Juzo requires that the garment be laundered (for hygienic reasons) prior to returning to the original place of purchase.

Fabric

Outer Shell - 88% Nylon, 12% Spandex / Inner Shell - 100% Neoprene

Instructions for Use

It is recommended that you consult with your medical professional and fitter for proper guidelines on the use and application of this product. If you experience any discomfort, immediately remove the short-stretch compression wrap and consult with your medical professional.

- Juzo Compression Wrap is designed to be reversible. The half circle hook fasteners need to be attached to the garment on the opposing side of the color to be worn.
- Juzo Compression Cotton Stockinette option provides a sanitary layer between the garment and the skin and/or any wound dressing.

Indications

Primary & secondary varices, chronic venous insufficiency, post-thrombotic syndrome, venous stasis ulcers, pre and post-vein procedure and lymphatic edema.

Contraindications

Absolute: peripheral arterial disease, decompensated heart diseases, septic phlebitis, phelgmasia coerulea dolens.

Relative: Leg infection, intolerance to compression or garment fabric, sensory disturbance to the limb.

Cautions

No side effects are currently known. However, if a negative result should occur (i.e. skin irritations) from the use of our products, please discontinue use of the garment and consult with a medical professional. The manufacturer is not liable for damage /injuries which have been caused by improper handling or misuse.

Steps to Putting on the Foot Wrap

1. The Juzo Compression Wrap comes with a compression cotton stockinette that can be worn optionally.

If using the Juzo Compression Cotton Stockinette, place the garment covering the area from the metatarsals (toes) to the ankle, with the cuffed edge closest to the metatarsals (toes)

2. Attach the hook fasteners to the straps on the side opposite of the color to be worn. For example, if beige is the desired color, then the fasteners should be applied to the black side of the garment.

2

Steps to Putting on the Foot Wrap

3. Gently fold the fasteners over onto themselves to prevent the fasteners from inadvertently attaching to the garment. Each fastener can then be unfolded before attaching the strap to the garment.

3

Steps to Putting on the Foot Wrap

4. Slide the wrap on your foot so that the heel fits into the heel cup. The seam at the back of the garment, above the heel opening should be centered at the back of the heel. If being worn with the Juzo Compression Wrap Calf Segment, it should be applied second and slightly overlap the bottom part of the calf segment.

4

Steps to Putting on the Foot Wrap

5. Simultaneously pull the bottom strap, nearest the toes, and the bottom EZ Pull until the strap is comfortably snug.

5

Steps to Putting on the Foot Wrap

6. After each strap is initially fastened, go back and pull the second and third straps again until they are both comfortably tightened around the heel and instep area.
7. Pull the top strap, nearest the ankle, and the top EZ Pull until it is comfortably tightened around the ankle.

Steps to Putting on the Foot Wrap

8. The wrap should have a snug fit but not painfully tight. If you feel any pain or discomfort or the wrap feels too loose, then repeat steps 6 through 8, adjusting the level of pressure that is applied

Note: If you are using a Juzo Compression Wrap Calf Segment, see care and application directions for use.

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®Juzo, FiberSoft & Keep Moving are registered trademarks of Julius Zorn, Inc. © 2018 Julius Zorn Inc.
Julius Zorn, Inc. • 3690 Zorn Dr.
JUZO CANADA LTD.