Visit [www.juzousa.com](http://www.juzousa.com) or [www.juzo.ca](http://www.juzo.ca) for more information.

Juzo USA  
P.O. Box 1088  
Cuyahoga Falls, OH 44223  
1 888-255-1300  
[www.juzousa.com](http://www.juzousa.com)  
[support@juzousa.com](mailto:support@juzousa.com)

Juzo Canada Ltd.  
1100 Burloak Dr. Ste. 300  
Burlington, ON L7L 6B2  
1 888-255-1300  
[www.juzo.ca](http://www.juzo.ca)  
[support@juzo.ca](mailto:support@juzo.ca)

® Juzo, Slippie Gator, Keep Moving, Fiber Soft & Signature Collection are registered trademarks of Julius Zorn, Inc. Lycra is a registered trademark of Invista Inc. Tyvek is a registered trademark of Dupont. © 2021 Julius Zorn, Inc.

Freedom in Motion

## Care & Application

### Compression Stockings & Pantyhose

### Instructions for Use

It is recommended that you consult with your medical professional and fitter for proper guidelines on the use and application of this product. As therapy progresses or new symptoms develop, it may be necessary to take new measurements for a new compression stocking.

It is also recommended that you put your compression stockings on first thing in the morning every day. Rough fingernails, jagged rings, etc. can damage fabric; use a Juzo latex free donning glove to protect garment.

### Indications

Primary & secondary varices, pregnancy varices and or edema, venous & chronic venous insufficiency, post superficial thrombophlebitis, prophylaxis of thrombosis & the post-thrombotic conditions, post-thrombotic syndrome, post-operative & post-traumatic edema, lymphatic edema & lipedema, elephantiasis, angiodyplasia, pre & post vein procedure, arthrogenous congestive syndrome, paresis & partial paresis of extremities, prophylaxis of venous disorders, heavy & fatigued legs, atrophy blanche, dermatosclerosis, post ulcer, severe tendency to edema & treatment of burn scars.

### Contraindications\*

Absolute: peripheral arterial disease, decompensated heart diseases, septic phlebitis, phlegmasia coerulea dolens.  
Relative: concomitant dermatoses, intolerance to compression stocking fabric, sensory disturbance of the limb, advanced peripheral neuropathy, primary chronic arthritis.

### Cautions

**Juzo Silver Garments should be removed prior to MRI Scans & Cardiac Defibrillations. In the event any side effects occur (i.e. skin irritation, etc.) from the use of a garment, you should consult with a medical professional.**

**The manufacturer is not liable for damage or injury caused by improper handling or misuse of its garments.**

\*Contraindications according to the German Phlebology Society.

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## Freedom to Make the Most of Every Day

Riding a bike. Working in the garden. Taking a walk. Going to work. These are the keep moving benefits of Juzo compression garments, designed to improve your quality of life, not inhibit it. Family-owned since 1912, Juzo has led the compression therapy garment industry by creating high-quality products that offer you the freedom to maintain a more active lifestyle. Advancements like our FiberSoft technology, seasonal colors, attractive styles and fun print designs ensure you move in comfort without compromising fashion. Thank you for choosing Juzo.

## Care and Maintenance

To maintain your Juzo stockings, we recommend the following:

- Hand or machine wash warm daily to maximize therapeutic effectiveness and prolong the life of your Juzo garment.
- When machine washing and drying, use the lowest permanent press settings. A garment bag is strongly recommended.
- Wash dark colors separately from other clothing.
- Use Juzo washing detergent or another mild soap.
- DO NOT USE: bleach, chlorine, fabric softeners, or other laundry additives. These may damage your garment.
- Please refer to the washing symbols on the label in your Juzo stocking.
- Store your stockings in a dry place at room temperature.
- Largely resistant to sun, ointments, lotions, and perspiration.

## Return Policy

Should it become necessary to return your garment for evaluation, for hygienic reasons Juzo requires that the garment be laundered prior to returning to the original place of purchase.

## Fabric

Details are provided on label sewn inside your compression stocking.

## Product Registration

Please take a few minutes to register your product online to enable us to provide you with best customer service possible, confirm ownership, help us to continue to improve product development and notify you of important product updates.

Register your Juzo product online at:

[www.juzoproducts.com](http://www.juzoproducts.com)

## Useful tips to help with your compression

### stocking: For knee-highs, thigh-highs and pantyhose:

- 1 Turn the upper half of the stocking inside out to the heel by gently sliding the stocking down over the foot portion.
- 2 Slip your foot into the foot portion.
- 3 Gently ease the remaining part of the stocking over the heel.
- 4 Use the palms of your hands to massage the stocking upwards towards the knee. NOTE: Never pull at the upper edge of the stocking.
- 5 Follow the same procedure for the opposite leg, taking care that the stitch rows always run vertically.
- 6 For thigh-highs and pantyhose, position both stockings, as described above, until just below the knee. NOTE: If putting on a thigh-high stocking with silicone border, never pull the stocking up from the silicone border.
- 7 Gently ease the stocking onto the thigh. Work the stocking upwards into place using the palms of your hands until the stocking reaches the groin area. NOTE: If putting on a thigh-high stocking with silicone border, never pull the stocking up from the silicone border.
- 8 For pantyhose, ease the body part of the pantyhose to your waist. The seams should run vertically up the front of the garment.

Always use Juzo latex free donning gloves when applying compression garments.

## Enclosed with your JUZO OPEN TOE compression stocking is a complimentary Slippie designed to help you put on your open toe compression garment.

Follow the below instructions.

- 1 Insert half the Juzo Slippie into the top of the compression stocking.
- 2 Place the compression stocking w/ the inserted Juzo Slippie lengthwise on the Easy Pad. (Easy Pad is not included)
- 3 Slip your foot into the Juzo Slippie while it is inside the compression stocking.
- 4 Put weight onto the entire foot. Slide it forward and lift your toes slightly during the movement. The compression stocking will slide over the heel.
- 5 Repeat the last step until your entire foot is in the stocking.
- 6 When the heel of the stocking is in place adjust the foot portion by sliding your foot sideways to the left or right on the Juzo Easy Pad.
- 7 Place the heel of your opposite foot on the tip of the Juzo Slippie. Pull our leg back removing Slippie from inside the stocking.
- 8 While wearing Juzo donning gloves, work the stocking further up the leg. Constantly correct the fit of the compression stocking and evenly distribute the material.

For more information about Juzo donning aids, go to: [www.juzousa.com](http://www.juzousa.com) • [www.juzo.ca](http://www.juzo.ca)

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